



Laboratory And Biological Services

M • I • C • R • O / A • N • A • L • Y • T • I • C • A • L C • O • N • S • U • L • T • A • N • T • S

LABORATORY SAMPLE REPORT: CERTIFICATE OF ANALYSIS

NANOCARE SA
600 CHRIS HANI RD
GREENWOOD PARK, DURBAN NORTH
UNIT 5
4051
082 495 3279
082 775 2282

TEST REP NO : 015/22
ISSUE NO : 001
REPORT DATE : 11/08/22
DATE REC : 11/07/22
DATE TESTED : 17/07/22
PAGE : 1 OF 2

ATTENTION: SHAUN VAN WYNGAARD

1. CONDITION OF SAMPLE ON RECEIPT PRIOR TO TESTING: ACCEPTABLE

2. DESCRIPTION OF SAMPLE: NANOCHLOR/NANOPURE-3.25G DISINFECTANT TABLET 52% NaDCC-0.008%
QUANTUM NANO SILVER 1x3.25G=1000ppm/1LT WATER
BN:18901 DATE:01/05/2022 EXPIRY:02/05/2025
PRODUCT BRAND NAME : NANOCARE

3. TEST REQUESTED AND OBJECTIVE OF THE TESTING:

Swab testing of treated surfaces, to determine the residual bactericidal activity of the sample/product at a 1,7 and 14 day period.

4. TEST METHOD REFERENCED:

SANS 5763: Testing of environmental swabs for the T.A.C. (Total Aerobic Count of the recoverable organisms present).

5. TEST PROCEDURE:

-A melamine surface, cleaned and disinfected with 70% ethanol, was used for the testing.
Ethanol was given sufficient time to dry and evaporate before any testing commenced.

-Surfaces were also swabbed for control testing purposes and to ensure levels of sterility during testing.

-All surfaces were marked for identification, with sufficient spacing to ensure completion of all tests.

-The surface was then sanitised/disinfected with the product (Nanochlor/Nanopure disinfectant tablet) in a solution of tap water, the test area was disinfected by the spray and wipe method.

-Surface area of 10cm² was then inoculated with the E.coli organism, using 1ml of the inoculum preparation (organism count used per ml $\pm 1.0 \times 10^9$ cfu/ml).

-After organism exposure and contact time of 3 hours, the surfaces were then swabbed for T.A.C, recovery of organisms, still present.

-The above test was carried out at the 1, 7 and 14 day period.

NOTE: LABS VALUES OUR CLIENT'S COMMENTS. PLEASE DO NOT HESITATE TO OFFER SUGGESTIONS OR COMMENTS THAT COULD IMPROVE OUR SERVICE

50 Hennie Alberts Street Brackenhurst 1448

P.O. Box 145406 Brackengardens 1452 • Tel: 011 867-2059/0593 • Fax: 011 867-3281 • E-mail labserve@mweb.co.za

Reg. No. CK 91/08374/23 • Members: **R.G. Khoury, **C.A. Khoury,

**Approved signatories: **E. Mohlaping

SEE OVERLEAF FOR TERMS AND CONDITIONS

cfu/g/ml: Colony forming units per gram or per ml of sample tested

*Not an accredited method

ND: Not detected



Laboratory And Biological Services

M • I • C • R • O / A • N • A • L • Y • T • I • C • A • L C • O • N • S • U • L • T • A • N • T • S

LABORATORY SAMPLE REPORT: CERTIFICATE OF ANALYSIS

TEST REP NO : 015/22

PAGE : 2 OF 2

6. TEST RESULTS: SPRAY AND WIPE METHOD, 3 HOUR EXPOSURE TIME

Organism: E.coli ATCC 8739

DAY 1

Viable Organisms Added 0 Time	Control Count	Spray and wipe Surface	% Kill Rate
6.0×10^8	5.0×10^6	4.0×10^4	99,20

DAY 7

Viable Organisms Added 0 Time	Control Count	Spray and wipe Surface	% Kill Rate
4.0×10^9	1.0×10^8	3.0×10^1	99,99

DAY 14

Viable Organisms Added 0 Time	Control Count	Spray and wipe Surface	% Kill Rate
3.0×10^9	9.0×10^6	0.0	100

7. INTERPRETATION OF RESULTS

The actives in the product (Nanochlor/Nanopure disinfectant tablet) within the prescribed exposure time of 3 hours has reached a 100% kill rate of the E.coli organisms present, with a range of between 99.20 and 100%. This was done using a spray and wipe method of application to the melamine surface at the 1, 7 and 14 day period.

COMMENTS:

REPORTED BY: R. KHOURY
MICROBIOLOGIST

REPORTED BY: E. MOHLAPING
LAB MANAGER

The above results reflect the analysis only of the sample in its conditions received. This report may not be reproduced except in full and without the written permission of LABS. This document is issued pursuant to service(s) provided within the scope of the clients instructions subject to the standard Terms and Conditions of Business of Laboratory and Biological Services CC(a copy of which is available on request) except that liability shall not exceed a sum equal to the fee paid for the said service and claims in respect to this document will be considered only if based upon proven gross negligence. This is not intended to relieve the parties to any sales contract from their contractual obligations.

NOTE: LABS VALUES OUR CLIENT'S COMMENTS. PLEASE DO NOT HESITATE TO OFFER SUGGESTIONS OR COMMENTS THAT COULD IMPROVE OUR SERVICE

50 Hennie Alberts Street Brackenhurst 1448

P.O. Box 145406 Brackengardens 1452 • Tel: 011 867-2059/0593 • Fax: 011 867-3281 • E-mail labserve@mwweb.co.za

Reg. No. CK 91/08374/23 • Members: **R.G. Khoury, **C.A. Khoury,

**Approved signatories: **E. Mohlaping

SEE OVERLEAF FOR TERMS AND CONDITIONS

cfu/g/ml: Colony forming units per gram or per ml of sample tested

*Not an accredited method

ND: Not detected



Laboratory And Biological Services

M • I • C • R • O / A • N • A • L • Y • T • I • C • A • L C • O • N • S • U • L • T • A • N • T • S

LABORATORY SAMPLE REPORT: CERTIFICATE OF ANALYSIS

NANOCARE SA
600 CHRIS HANI RD
GREENWOOD PARK, DURBAN NORTH
UNIT 5
4051
082 495 3279
082 775 2282

TEST REP NO : 016/22
ISSUE NO : 001
REPORT DATE : 11/08/22
DATE REC : 11/07/22
DATE TESTED : 17/07/22
PAGE : 1 OF 2

ATTENTION: SHAUN VAN WYNGAARD

1. CONDITION OF SAMPLE ON RECEIPT PRIOR TO TESTING: ACCEPTABLE

2. DESCRIPTION OF SAMPLE: NANOCHLOR/NANOPURE-3.25G DISINFECTANT TABLET 52% NaDCC-0.008%
QUANTUM NANO SILVER 1x3.25G=1000ppm/1LT WATER
BN:18901 DATE:01/05/2022 EXPIRY:02/05/2025
PRODUCT BRAND NAME : NANOCARE

3. TEST REQUESTED AND OBJECTIVE OF THE TESTING:

Swab testing of treated surfaces, to determine the residual bactericidal activity of the sample/product at a 1,7 and 14 day period.

4. TEST METHOD REFERENCED:

SANS 5763: Testing of environmental swabs for the T.A.C. (Total Aerobic Count of the recoverable organisms present).

5. TEST PROCEDURE:

-A melamine surface, cleaned and disinfected with 70% ethanol, was used for the testing.
Ethanol was given sufficient time to dry and evaporate before any testing commenced.

-Surfaces were also swabbed for control testing purposes and to ensure levels of sterility during testing.

-All surfaces were marked for identification, with sufficient spacing to ensure completion of all tests.

-The surface was then sanitised/disinfected with the product (Nanochlor/Nanopure disinfectant tablet) in a solution of tap water, the test area was disinfected by the spray and wipe method.

-Surface area of 10cm² was then inoculated with the Staph Aureus organism, using 1ml of the inoculum preparation (organism count used per ml $\pm 1.0 \times 10^9$ cfu/ml).

-After organism exposure and contact time of 3 hours, the surfaces were then swabbed for T.A.C, recovery of organisms, still present.

-The above test was carried out at the 1, 7 and 14 day period.

NOTE: LABS VALUES OUR CLIENT'S COMMENTS. PLEASE DO NOT HESITATE TO OFFER SUGGESTIONS OR COMMENTS THAT COULD IMPROVE OUR SERVICE

50 Hennie Alberts Street Brackenhurst 1448

P.O. Box 145406 Brackengardens 1452 • Tel: 011 867-2059/0593 • Fax: 011 867-3281 • E-mail labserve@mweb.co.za

Reg. No. CK 91/08374/23 • Members: **R.G. Khoury, **C.A. Khoury,

**Approved signatories: **E. Mohlaping

SEE OVERLEAF FOR TERMS AND CONDITIONS

cfu/g/ml: Colony forming units per gram or per ml of sample tested

*Not an accredited method

ND: Not detected



Laboratory And Biological Services

M•I•C•R•O / A•N•A•L•Y•T•I•C•A•L C•O•N•S•U•L•T•A•N•T•S

LABORATORY SAMPLE REPORT: CERTIFICATE OF ANALYSIS

TEST REP NO : 016/22
PAGE : 2 OF 2

6. TEST RESULTS: SPRAY AND WIPE METHOD, 3 HOUR EXPOSURE TIME

Organism: Staph Aureus ATCC 6538

DAY 1

Viable Organisms Added 0 Time	Control Count	Spray and wipe Surface	% Kill Rate
2.0×10^8	4.0×10^7	2.3×10^7	42,50

DAY 7

Viable Organisms Added 0 Time	Control Count	Spray and wipe Surface	% Kill Rate
7.0×10^9	9.0×10^8	7.0×10^4	99,99

DAY 14

Viable Organisms Added 0 Time	Control Count	Spray and wipe Surface	% Kill Rate
7.3×10^9	2.6×10^9	3.9×10^3	99,99


7. INTERPRETATION OF RESULTS

The actives in the product (Nanochlor/Nanopure disinfectant tablet) within the prescribed exposure time of 3 hours has reached a 99,99% kill rate of the Staph Aureus organisms present, with a range of between 42,50% and 99,99%. This was done using a spray and wipe method of application to the melamine surface at the 1, 7 and 14 day period.

COMMENTS:



REPORTED BY: R. KHOURY
MICROBIOLOGIST



REPORTED BY: E. MOHLAPING
LAB MANAGER

The above results reflect the analysis only of the sample in its conditions received. This report may not be reproduced except in full and without the written permission of LABS. This document is issued pursuant to service(s) provided within the scope of the clients instructions subject to the standard Terms and Conditions of Business of Laboratory and Biological Services CC(a copy of which is available on request) except that liability shall not exceed a sum equal to the fee paid for the said service and claims in respect to this document will be considered only if based upon proven gross negligence. This is not intended to relieve the parties to any sales contract from their contractual obligations.

NOTE: LABS VALUES OUR CLIENT'S COMMENTS. PLEASE DO NOT HESITATE TO OFFER SUGGESTIONS OR COMMENTS THAT COULD IMPROVE OUR SERVICE

50 Hennie Alberts Street Brackenhurst 1448

P.O. Box 145406 Brackengardens 1452 • Tel: 011 867-2059/0593 • Fax: 011 867-3281 • E-mail labserve@mwweb.co.za

Reg. No. CK 91/08374/23 • Members: **R.G. Khoury, **C.A. Khoury,

**Approved signatories: **E. Mohlaping

SEE OVERLEAF FOR TERMS AND CONDITIONS

cfu/g/ml: Colony forming units per gram or per ml of sample tested

*Not an accredited method

ND: Not detected